



Masks and Prongs

[Medijet® · Miniflow®]

Masks



REF see table >>

The masks have to be used in combination with the Medijet® or the Miniflow®. All masks are made of silicone. For choosing the correct size use the medin® measuring tape.

technical data & facts

Size	x-small	small	medium	large	x-large
REF	1200-08	1200-04	1200-05	1200-06	1200-07
Mask width nose [mm]	14	18,5	21,0	23,0	26,0
Packaging unit	10	10	10	10	5

Prongs



REF see table >>

The prongs have to be used in combination with the Medijet® or the Miniflow®. All prongs are made of silicone. For choosing the correct size use the medin® measuring tape.

technical data & facts

Size	x-small	small	medium	large	x-large	mediumwide	large wide
REF	1200-01	1200-21	1200-02	1200-22	1200-03	1200-32	1200-33
Nostril Ø [mm]	3	3,5	4,1	4,75	5,5	3,7	5
Packaging unit	10	10	10	10	10	10	10